



Drawing





Material 01 Body - PE



Layout/Inner box

Inner box size: 145x60x110 mm
Inner box qty: 100pcs



Layout/Soft blister

Blister size:
54x26,6 mm


Printing size:
39x20 mm

Printing color:
PANTONE 200/288

Print without frame



Layout/Sticker on carton



Printing size:
180x140 mm

Printing color:
PANTONE 200/288

2 stickers glued on
shipping carton

For Combistopper BLUE is used blue colour and relevant texts

Packaging

corrugated cardboard, without iron pins, no plastic strips

Size: 455x365x360 mm

Volume: 0,05978 cbm

Ctn: 5000 pcs

Palette: 100.000 pcs (20 ctn)

Weight ctn: 8 kg

1.2 x 0.8 x 2.0m

Sterilization

EO sterilization

Guaranty

5 years

PRODUCTION LIST

CPH-PrL-CMS rev.1/02-2022